

Nucu^{Syke} manual

Table of Contents

1)	Contents of the package	3
2)	How to use Nucu	4
3)	How to clean Nucu	5
4)	Assembling the Nucu after cleaning	7
5)	Technical specifications	7
6)	Support	7

1) CONTENTS OF THE PACKAGE

Confirm that the following parts are included:



Nucu



Nucu base plate (inside the Nucu)



AC adapter



USB - Micro USB Cable

Note:

1. If any part is damaged, do not use it. Contact Nucu customer service:
 - info@nucu.fi
2. In hospital use check that the AC adapter has the “Approved Medical Device” notation



2) HOW TO USE NUCU

Place the Nucu on a flat surface and put the baby in the nest. Observe how the baby reacts to Nucu. Do not leave the baby unattended on the Nucu.

Functions	What to do
Power On	Plug the AC adapter to the wall outlet
Power Off	Unplug the AC adapter from the wall outlet

3) HOW TO CLEAN NUCU

Nucu's fabric parts are machine washable. Base plate can be cleaned with a disinfectant. Wash Nucu when necessary.

Fabric material: 100% cotton. Care symbols:     

Steps:

1) Open the

a) padding's zipper

b) bottom section zipper



2) Pull out the

a) base plate and mattress

b) padding



3) Wipe the base plate with a disinfectant



4) Wash the fabric parts in the washing machine



5) Allow the fabric to dry well

4) ASSEMBLING THE NUCU AFTER CLEANING

You can assemble the Nucu after cleaning when all parts are dry.

Steps:

- 1) Place the mattress inside the base section
- 2) Slide the base plate under the mattress
- 3) Place the padding in its place
- 4) Close the base section's and padding's zippers carefully

5) TECHNICAL SPECIFICATIONS

- Recommended operating temperature -10°C – 40°C
- Protection against dust ingress and splash water. IP54
- Nucu is designed and manufactured in Finland



6) SUPPORT

- info@nucu.fi
- www.nucu.fi/en/support